REF 925702

08.22

de

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NANOCONTROL NANOTURB

Methode:

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08.22

en

NANOCOLOR® Analytical Quality Control

NANOCONTROL NANOTURB

Method:

Turbidity standards for nephelometric turbidity measurements (primary standards for verification and calibration in accordance to ISO 7027)

Suitable for **NANOCOLOR®** spectrophotometers (*UV/VIS*, *VIS*, *UV/VIS II*, *VIS II* and *Advance*) and for compact photometer *PF-12^{Plus}*

Range:	1–400 NTU
Wavelength:	860 nm
Storage temperature:	20–25 °C

Contents of reagent set:

1 tube with test solution 1 NTU	1 tube with test solution 100 NTU
1 tube with test solution 4 NTU	1 tube with test solution 400 NTU
1 tube with blank value "NULL"	1 lot certificate

Note:

The test solutions are primary standards (NIST traceable standards).

Hazard warning:

Information regarding safety can be found on the box' label and in the safety data sheet. You can download the SDS from www.mn-net.com/SDS.

Procedure:

Measuring of turbidity standards:

The turbidity standards can be measured in the instruments *UV/VIS*, *VIS*, *UV/VIS II*, *VIS II* and *PF-12^{Plus}*. For measurement of the turbidity method 906 can be used at 860 nm. A nephelometric measurement is performed at a 90° angle. The result is displayed in NTU.

Instrument calibration:

Procedure with **NANOCOLOR** *UV/VIS* and *VIS*:

1. Choose **Turbidity calibration** through the menu **Calibration**.
2. Enter the number of calibration points and confirm with **OK**.
3. Enter the standard values and confirm with **OK**.
4. By pressing **Measure** the calibration is started.
5. Perform a zero measurement with distilled water (blank value), then run a turbidity calibration according to the instructions in the display using the turbidity standards. The measurements are completed automatically as soon as the respective cuvette is inserted.
6. Confirm the calibration with **OK**.

Procedure with **NANOCOLOR** *UV/VIS II* and *VIS II*:

1. Choose **Turbidity calibration** through the main menu icon and the menu **Calibration**.
2. Choose the option „NANOTURB, REF 925702 (> 1 NTU)“, the calibration will be started by pressing the **Measure** icon.
3. Perform a zero measurement with distilled water (blank value), then run a turbidity calibration according to the instructions in the display using the turbidity standards. The measurements are completed automatically as soon as the respective cuvette is inserted.
4. The storage comes automatically after performing the last measurement.

Procedure with **PF-12^{Plus}**:

1. Choose directly method 906 or use the menu to select **Basic functions**.
2. By pressing the button **Null / Zero** the turbidity calibration is started.
3. Perform a zero measurement with distilled water (blank value), then run a turbidity calibration according to the instructions in the display using the turbidity standards. Press **M** to perform the measurements.
4. The storage comes automatically after performing the last measurement.

Procedure with **NANOCOLOR® Advance**:

1. Choose **Turbidity calibration** through the menu **Calibration** from the third page of the user interface.
2. Choose the option „14 mm Küvette (NANOTURB, REF 925702)“. The calibration will be started by pressing the **Measure** icon.
3. Perform a zero measurement with distilled water (0 NTU), then run a turbidity calibration according to the instructions in the display using the turbidity standards. The measurements are performed automatically as soon as the respective cuvette is inserted.
4. The storage comes automatically after performing the last measurement.

Note:

Tubes must be clean. If necessary, clean them with a dry and clean cloth. Do not use scratched tubes. For turbidity calibrations we recommend to use the cuvette slot cover.

Disposal:

Information regarding disposal can be found in the safety data sheet. You can download the SDS from www.mn-net.com/SDS.

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