# REF 918978 03.24 NANOCOLOR® NanOx Metall

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The following determinations can be performed with the NANOCOLOR® system after oxidation with NanOx Metal:

Test 0-98 Aluminium 07*	Test 1-36 Iron
Test 0-2 Aluminium*	Test 0-37 Iron 3
Test 0-14 Cadmium 2	Test 0-71 Nickel 4
Test 0-24 Chromate 5 (Chromium)	Test 1-62 Nickel
Test 1-25 Chromate (Chromium)	Test 0-76 Phosphate 1
Test 1-51 Cobalt	Test 0-81 Phosphate 5
Test 1-53 Copper	Test 0-80 Phosphate 15
Test 0-54 Copper 7	Test 0-55 Phosphate 45
Test 0-53 Copper 5	Test 0-79 Phosphate 50
	Test 1-95 Zinc

Test 0-96 Zinc 4

Test 0-42 Zinc 6

\*only for microwave decomposition

### Nitrogen

For determination of total nitrogen NANOCOLOR® NanOx N (REF 918979) is recommended.

### General notes for the decomposition procedure:

Please use the orange measuring spoon only for the decomposition reagent. Don't use the neutralization reagent for following tests: chromium test 0-24/1-25, iron test 1-36/0-37, copper test 1-53/0-54 and phosphorus test 0-79. When performing copper tests 1-53 and 0-54 use twice as much reagent R1 (test 1-53) or R2 (test 0-54) respectively (REF 918903). Dilute the sample until the measured value is within the measuring range. For waters of unknown concentrations we recommend that you perform the test with very different dilutions (e.g. 1+9, 1+99) until the last dilution confirms the previous value. For samples which consume large amounts of oxidizing substances (e.g. for COD values above 1000 mg/L O<sub>2</sub>), decomposition can be incomplete, or the safety diaphragm in the pressure vessel can be destroyed. In such cases repeat the decomposition with a diluted sample solution. The power of commercial microwave ovens can fluctuate, thus it is necessary to check the microwave power from time to time and eventually optimize the microwave time.

### Decomposition in the heating block:

#### Decomposition

##### Requisite accessories:

NANOCOLOR® heating block

NANOCOLOR® test tubes 16 mm OD (REF 91680)

Pipet **6.0 mL** sample solution into an empty reaction tube 16 mm OD (REF 91680), add **1 level orange measuring spoon NanOx Metal decomposition reagent**, close and shake thoroughly. Place the reaction tube into the NANOCOLOR® heating block and heat 30 min at 120 °C or 1 h at 100 °C. Remove tube from the heating block, shake lightly and let cool. **The decomposition solution must be clear and colourless.** Otherwise another decomposition is necessary. Open the reaction tube and test the decomposition solution for peroxides using QUANTOFIX® Peroxide 25 test sticks (REF 91319). If peroxides are present, close tube and heat again without further addition of NanOx Metal decomposition reagent. Remove from the heating block and let cool for about 10 min. Turn the reaction tube on its head once, open it and again test for peroxides. Then carefully (evolution of gas) add **3 level white microspoons NanOx Metal neutralization reagent**, close and shake thoroughly. The pH value must be pH = 3–7, otherwise add more neutralization reagent. The contents can now be used as sample solution for the above-mentioned phosphorus and metal determinations. Please note that for standard tests the measured result has to be multiplied with the dilution factor depending upon the used amount of sample.

### Decomposition XL:

#### Requisite accessories:

NANOCOLOR® heating block (REF 919350.1)

NANOCOLOR® decomposition vessel 22 mm OD (REF 91622)

Pipet **17 mL** sample solution into an empty reaction tube 22 mm ID (REF 91622), add **3 level orange measuring spoon NanOx Metal decomposition reagent**, close and shake thoroughly. Place the reaction tube into the NANOCOLOR® heating block and heat 30 min at 120 °C or 1 h at 100 °C. Remove tube from the heating block, shake lightly and let cool. **The decomposition solution must be clear and colorless.** Otherwise another decomposition is necessary. Turn the reaction tube upside down once, open the tube and test the decomposition solution for peroxides using QUANTOFIX® Peroxide 25 test sticks (REF 91319). If peroxides are present, close tube and heat again without further addition of NanOx Metal decomposition reagent. Remove from the heating block and let cool for about 10 min. Turn the reaction tube upside down once, open it and again test for peroxides. If it is free of peroxides, then carefully (evolution of gas) add **3 level black spoons NanOx Metal neutralization reagent**, close and shake thoroughly. The pH value must be pH = 3–7, otherwise add more neutralization reagent. The contents can now be used as sample solution for the above-mentioned phosphorus and metal determinations. Please note that for standard tests the measured result has to be multiplied with the dilution factor depending upon the used amount of sample. Please remove residual solution from the decomposition vessel after use and clean it with a bottle brush and distilled water.

## &lt;h

<b>Decomposition:</b>	Pipet 60 $\mu$ L test sample and 5.94 mL dist. water into a 14 mm ID reaction tube
<b>Decomposition XL:</b>	Pipet 170 $\mu$ L test sample and 16.8 mL dist. water into the decomposition vessel

#### Alternative procedure for microwave decomposition:

##### Requisite accessories:

Microwave oven with digital display (time setting in seconds), PTFE pressure vessel 20 mL for decomposition and safety expansion reservoir (REF 91640)

Pipet 10 mL sample solution into the decomposition vessel, add 2 level orange measuring spoons Nanox Metal decomposition reagent, close and shake thoroughly. Place pressure vessel on the outer rim of the turntable in the microwave oven and heat 23 s at 900 VA or 28 s at 750 VA (always use the highest setting of the respective instrument). Remove vessel from the microwave oven and let cool for about 10 min. Turn the pressure vessel on its head once and then open it with caution. **The decomposition solution must be clear and colorless.** Otherwise another decomposition is necessary. Test the decomposition solution for peroxides using QUANTOFIX® Peroxide 25 test sticks (REF 91319). If peroxides are present, close vessel, and microwave again without further addition of Nanox Metal decomposition reagent. Remove from the microwave oven and let cool for about 10 min. Turn the pressure vessel on its head once, open it and again test for peroxides. Then carefully (evolution of gas) add 6 level white microspoons Nanox Metal neutralization reagent, close and shake thoroughly. The pH value must be pH = 3–7, otherwise add more neutralization reagent. The contents can now be used as sample solution for the above-mentioned phosphorus and metal determinations. Please note that for standard tests the measured result has to be multiplied by 2 due to application of half the sample volume.

##### Dilutions:

Sample 1+9 diluted	
<b>Microwave:</b>	Pipet 1.0 mL test sample and 9.0 mL dist. water into the decomposition vessel
Sample 1+99 diluted	
<b>Microwave:</b>	Pipet 0.1 mL test sample and 9.9 mL dist. water into the decomposition vessel

##### Interferences:

###### a) The safety diaphragm ruptures

- Material fatigue  
→ replace safety diaphragm and sealing disc
- Sample develops gases → dilute sample
- Delays in boiling  
→ use special Nanox boiling chips
- Wrong time setting → adjust time

###### b) Incomplete decomposition

- Decomposition incomplete due to high COD content  
Note: slightly nitrous odour when opening the decomposition vessel or **decomposition solution not clear**  
→ dilute into the measuring range, if necessary decompose several times
- other interferences, which cannot be eliminated by decomposition with Nanox Metal

##### Time adjustment for different microwave ovens:

For decomposition with Nanox Metal the microwave times have to be adjusted for different instruments. All instruments should be operated at the highest setting, and the time is adjusted to the second.

##### Example for calculation:

Microwave energy 900 VA → 24 s heating time  
Microwave energy 750 VA → X s heating time  
 $X = 900 / 750 \times 24 \text{ s} = 28.8 \text{ s} \rightarrow 29 \text{ s}$

Please note that the actual power of the microwave oven may deviate from the nominal value. In these cases test and adjust the time as described below.

##### Correction and optimization of microwave times:

The output of microwave ovens can fluctuate even between instruments of the same model or it can change with time. For the decomposition with Nanox Metal, however, an exact dosing of energy is required. In order to optimize the microwave time, always choose the highest setting of the respective instrument (max. 1000 VA), because this is not interrupted by time intervals. For a microwave oven with 900 VA output the microwave time will be about 24 s. The following procedure is a very reliable and simple method to determine the microwave energy:

Fill the included 100 mL Erlenmeyer flask with exactly 100 mL water (of ambient temperature) and measure the temperature with a thermometer (e.g. 21 °C). Now place the filled Erlenmeyer flask (without thermometer) on the outer rim of the microwave turntable (same position as the decomposition vessel) and heat for the calculated time. Immediately remove the Erlenmeyer flask and measure the water temperature after slight stirring (e.g. 53 °C). The temperature difference must be 30 °C. In our example the difference is 53 - 21 = 32 °C. The heating time has to be cor-

rected for this deviation:

$$\text{Nominal difference/actual difference} \times \text{time} = \text{corrected time}$$

$$30 / 32 \times 24 \text{ s} = 22.5 \text{ s} \rightarrow 23 \text{ s}$$

Check this time setting again following the above procedure.

##### General remarks:

Always place the decomposition vessel on the same position of the microwave turntable. For easier handling we recommend that you fill the safety expansion reservoir with filter flocs or pulp, covered with a filter circle, and remove the exhaust tubing from the safety expansion reservoir. After 5 decompositions the filter flocs (REF 481110) or pulp have to be replaced. If the safety diaphragm ruptures, always dry the safety reservoir prior to the next run, or replace the filter flocs or pulp. For sample solutions which tend to delays in boiling, we recommend using PTFE boiling chips (REF 91628). After decomposition remove the decomposition solution from the vessel. When the vessel is heavily soiled, or after continued use, clean the decomposition vessel with a bottle brush and dist. water. We recommend that from time to time you perform a blank decomposition with dist. water and Nanox Metal decomposition reagent, but without addition of neutralization reagent. As an alternative, you may fill the vessel with dist. water and some drops hydrogen peroxide 30% and let it stand for some hours.

##### REF

91312	QUANTOFIX® Peroxide 100
91319	QUANTOFIX® Peroxide 25
91640	PTFE pressure vessel 20 mL for decomposition with safety expansion reservoir, set of sealing discs and safety diaphragms
91643	PTFE pressure vessel 20 mL for decomposition
91640.1	PFA safety expansion reservoir
91646	PFA tubing with fitting
91644	pack of 20 PFA sealing discs
91645	pack of 10 PTFE safety diaphragms
91628	pack of 3 PTFE boiling chips
91680	pack of 20 reaction tubes
91636	rack for 15 reaction tubes and 2 decomposition vessels
91622	pack of 2 decomposition vessels 22 mm OD
481100	filter flocs, pack of 500 g
481110	filter flocs, pack of 1000 g

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