


NANOCOLOR® Sauerstoff 12

Oxygen 12/Oxygène 12/Oxígeno 12










Test 0-82
REF 985082

Methode / Method / Méthode / Método
0821 0.5–12.0 mg/L O₂


436 nm

Min	20–80 %		Max
0.5	2.8	9.7	12.0 mg/L



NANOCOLOR® Sauerstoff 12

Oxygen 12 / Oxygène 12 / Oxígeno 12



Test 0-82
REF 985082

Methode / Method / Méthode / Método
0821 0.5–12.0 mg/L O₂

436 nm

Min	20–80 %		Max
0.5	2.8	9.7	12.0 mg/L

