


NANOCOLOR® leicht freisetzbares Cyanid 04

easily liberated cyanide 04 / cyanure facilement libéré 04 / cianuro fácilmente liberado 04











Test 0-25

REF 985025

Methode / Method / Méthode / Método
0251 0.01–0.40 mg/L CN⁻


585 / 605 nm

Min	20–80 %		Max
0.01	0.09	0.32	0.40 mg/L CN ⁻



NANOCOLOR® leicht freisetzbares Cyanid 04

easily liberated cyanide 04 / cyanure facilement libéré 04 / cianuro fácilmente liberado 04



Test 0-25

REF 985025

Methode / Method / Méthode / Método

0251 0.01–0.40 mg/L CN⁻

585 / 605 nm

Min

0.01

20–80 %

0.09

0.32

Max

0.40 mg/L CN⁻

